GE Healthcare

# Gel Filtration Calibration Kits

Product booklet

Codes: 28-4038-41 Low Molecular Weight 28-4038-42 High Molecular Weight



## Page finder

1.	Legal	3	
2.	Handling	4	
	2.1. Safety warnings and precautions	4	
	2.2. Storage	4	
	2.3. Expiry	4	
3.	Components	5	
	3.1. Components in Gel Filtration Calibration kits		
	LMW and HMW	5	
	3.2. Common abbreviations	6	
4.	Introduction	7	
	4.1. Critical parameters	9	
5.	Procedure	10	
6.	5. Typical results		
7.	Important notes	20	
	7.1. The use of the calibration kits with denaturing		
	solvents	20	
	7.2. Dimer and oligomer formation in calibration		
	kit proteins	20	
	7.3. Electrophoresis calibration kits	21	
8.	Appendix	22	
9.	. Product information 2		

## 1. Legal

GE, imagination at work and GE monogram are trademarks of General Electric Company.

HiLoad, HiPrep, Sephacryl, Sephadex, Sepharose, Superdex, Superose, Tricorn, ÄKTA, and ÄKTAexplorer are trademarks of GE Healthcare companies.

**Tricorn Columns:** Tricorn column and components are protected by US design patents USD 500856, USD 506261, USD 500555, USD 495060 and their equivalents in other countries.

© 2006–2009 General Electric Company-All rights reserved. First published 2006.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

http://www.gelifesciences.com/protein-purification

GE Healthcare UK Limited. Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

## 2. Handling

# 2.1. Safety warnings and precautions

## Warning: For research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

All chemicals should be considered as potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing such as laboratory overalls, safety glasses and gloves. Care should be taken to avoid contact with skin or eyes. In the case of contact with skin or eyes, wash immediately with water. See material

safety data sheet(s) and/or safety statement(s) for specific advice.

### 2.2. Storage

The kit should be stored at 2°C to 8°C. It is recommended that each vial of protein is dissolved in a buffer with a pH of 6–8 and an ionic strength of  $\geq$  0.15 M (e.g 50 mM phosphate, 0.15 M NaCl, pH 7.2)

### 2.3. Expiry

For expiry details please see outer packaging.

## 3. Components

# 3.1. Components in Gel Filtration Calibration kits LMW and HMW

**Table 1.** The content of Gel Filtration Calibration Kit LMW (lowmolecular weight), code no. 28-4038-41.

Protein (weight per vial)	Molecular weight (M <sub>r</sub> )	Stoke´s Radius¹ (Å)	Source
Aprotinin (10 mg)	6500	NA	Bovine lung
Ribonuclease A (50 mg)	13 700	16.4	Bovine pancreas
Carbonic Anhydrase (15 mg)	29 000	NA	Bovine erythrocytes
Ovalbumin (50 mg)	44 000	30.5	Hen egg
Conalbumin (50 mg)	75 000	NA	Chicken egg white
Blue Dextran 2000			

 Table 2. The content of Gel Filtration Calibration Kit HMW (high molecular weight), code no. 28-4038-42.

Protein (weight per vial)	Molecular weight (M <sub>r</sub> )	Stoke´s Radius¹ (Å)	Source
Ovalbumin (50 mg)	44 000	30.5	Hen egg
Conalbumin (50 mg)	75 000	NA	Chicken egg white
Aldolase <sup>2</sup> (50 mg)	158 000	48.1	Rabbit muscle
Ferritin <sup>2</sup> (15 mg)	440 000	61.0	Horse spleen
Thyroglobulin (50 mg)	669 000	85.0	Bovine thyroid
Blue Dextran 2000			

- 1. Reference: CRC Practical Handbook of Biochemistry and Molecular Biology G.D. Fasman, ed., CRC Press, 1989, 601 pp.
- **2.** These proteins are supplied mixed with sucrose or mannitol to maintain stability and aid their solubility.

### 3.2. Common abbreviations

- **GF** = gel filtration
- GL = glass
- HMW = high molecular weight
- **HR** = high resolution
- K<sub>av</sub> = partition coefficient
- LMW = low molecular weight
- M<sub>r</sub> = relative molecular weight
- PC = precision column
- pg = prep grade
- **R**<sub>st</sub> = Stoke's radius of solute
- V<sub>0</sub> = void volume
- V<sub>c</sub> = geometric column volume
- V<sub>e</sub> = elution volume
- V<sub>t</sub> = total liquid volume

### 4. Introduction

Two **Gel Filtration Calibration Kits** are available for protein molecular weight determination by gel filtration. The **Low Molecular Weight Kit** contains 5 proteins with molecular weights in the range 6500 to 75 000 and Blue Dextran 2000. The **High Molecular Weight Kit** contains 5 proteins with molecular weights in the range 44 000 to 669 000 and Blue Dextran 2000.

The use of gel filtration chromatography for the determination of the molecular weight and size of proteins is well documented. The technique is based on the well-established ability of gel filtration media, such as Superdex<sup>™</sup>, Superose<sup>™</sup>, Sephacryl<sup>™</sup>, Sephadex<sup>™</sup> and Sepharose<sup>™</sup> to separate molecules according to size. Prepacked columns are available (Table 3) and can be run on chromatography systems such as, ÄKTA<sup>™</sup> design.

Molecular weight determination by gel filtration are carried out by comparing an elution volume parameter, such as  $K_{av}$  of the protein of interest, with the values obtained for several known calibration standards. In practice it is found that for homologous series of compounds a sigmoid relationship exists between their various elution volume parameters and the logarithm of their molecular weights.

The molecular weight of an unknown protein can be determined from the calibration curve (plot of  $K_{av}$  versus the logarithm of molecuar weight) once its  $K_{av}$  value is calculated from its measured elution volume.

For accurate determination of molecular weight, the calibration standards must have the same relationship between molecular weight and molecular size as the substance of interest. GE Healthcare Calibration Kits provide highly purified, wellcharacterized, globular protein standards for protein molecular weight determination.

Recommended Fractionation range sample volume (M, of globular (ul) durina Product  $V_c$  (ml) calibration protein) Superdex – Prepacked Tricorn™ and PC columns  $3 \times 10^3$  to  $7 \times 10^4$ Superdex 75 10/300 GL 24 100 Superdex 75 5/150  $3 \times 10^3$  to  $7 \times 10^4$ 3 12.5  $3 \times 10^3$  to  $7 \times 10^4$ Superdex 75 PC 3.2/30 24 10 Superdex 200 10/300 GL  $1 \times 10^4$  to  $6 \times 10^5$ 100 24  $1 \times 10^4$  to  $6 \times 10^5$ 3 Superdex 200 5/150 125 Superdex 200 PC 3.2/30  $1 \times 10^4$  to  $6 \times 10^5$ 2.4 10 Superdex prep grade – HiLoad™ prepacked columns HiLoad 16/60 Superdex 75 pg  $3 \times 10^3$  to  $7 \times 10^4$ 120 500 HiLoad 26/60 Superdex 75 pg  $3 \times 10^3$  to  $7 \times 10^4$ 320 1000 HiLoad 16/60 Superdex 200 pa  $1 \times 10^4$  to  $6 \times 10^5$ 120 500  $1 \times 10^4$  to  $6 \times 10^5$ HiLoad 26/60 Superdex 200 pg 320 1000 Superose – Prepacked Tricorn and PC columns  $1 \times 10^3$  to  $3 \times 10^5$ Superose 12 10/300 GL 24 100 Superose 12 PC 3.2/30  $1 \times 10^3$  to  $3 \times 10^5$ 24 10  $5 \times 10^3$  to  $5 \times 10^6$ Superose 6 10/300 GL 24 100  $5 \times 10^3$  to  $5 \times 10^6$ Superose 6 PC 3.2/30 2.4 10 Sephacryl – HiPrep™ prepacked columns HiPrep 16/60 Sephacryl S-100 HR  $1 \times 10^{3}$  to  $1 \times 10^{5}$ 120 500 HiPrep 26/60 Sephacryl S-100 HR  $1 \times 10^{3}$  to  $1 \times 10^{5}$ 320 1000 HiPrep 16/60 Sephacryl S-200 HR  $5 \times 10^3$  to  $2.5 \times 10^5$ 120 500  $5 \times 10^3$  to  $2.5 \times 10^5$ HiPrep 26/60 Sephacryl S-200 HR 320 1000 HiPrep 16/60 Sephacryl S-300 HR  $1 \times 10^4$  to  $1.5 \times 10^6$ 120 500 HiPrep 26/60 Sephacryl S-300 HR  $1 \times 10^4$  to  $1.5 \times 10^6$ 320 1000

 Table 3. Prepacked columns with differing media and separation ranges.

### 4.1. Critical parameters

- Select buffer with pH 6 to 8 and ionic strength  $\geq$  0.15 M (e.g. 0.15 M NaCl)
- Use prepacked columns
- Use recommended flow rate
- Use a sample volume that is 0.1 to 2% of total column volume
- Select a gel filtration medium with suitable fractionation range for sample (Table 3).
- Minimize the delay volume in the chromatography system from injection valve to column and from column to detector.

### 5. Procedure

#### 1) Selection of gel filtration column

Select a column with a fractionation range so that the expected molecular weight of your sample falls approximately in the middle of the range for that column. Superdex 200 can be used for a preliminary, quick approximation of the sample's molecular weight. A bed length of 30 to 60 cm is sufficient for most determinations.

#### 2) Equilibration of the column

If the column has been stored in 20% ethanol, wash the column first with 2 column volumes of distilled water. Equilibrate the column with 2 column volumes of buffer. A buffer with a pH of 6–8 and an ionic strength  $\geq$  0.15 M is suggested. A typical buffer is 50 mM phosphate in 0.15 M NaCl, pH 7.2.

#### 3) Choice of calibration kit proteins

Include Calibration Kit proteins of a higher molecular weight and of a lower molecular weight than that of the sample. The proteins listed in Table 4 may be mixed to give resolved peaks.

For calibration of	Mix a	Mix b	Mix c
Superdex 200	F+C+CA+R	Ald+O+R+Apr	
Superdex 75	C+CA+R+Apr	O+R+Apr	
Superose 6	T+Ald+CA+Apr	F+O+R	
Superose 12	F+O+R	Ald+CA+Apr	
Sephacryl 300	T+Ald+CA	F+O+R+Apr	
Sephacryl 200	Ald+CA+R+Apr	C+CA+R	O+CA+R
Sephacryl 100	Ald+CA+R+Apr	C+CA+R	O+CA+R

Table 4. Suitable protein mixtures for resolved peaks.

Apr - Aprotinin, R - Ribonuclease A, CA - Carbonic Anhydrase,

O - Ovalbumin, C - Conalbumin, Ald - Aldolase, F - Ferritin,

T - Thyroglobulin

#### 4) Individual protein preparation

We recommend that the proteins are dissolved in high concentration (20 mg/ml) and diluted with buffer before use. Dissolve the content of the vial in a buffer with a pH of 6–8 and an ionic strength of  $\geq$  0.15 M (e.g 50 mM phosphate, 0.15 M NaCl, pH 7.2). Ferritin and aldolase are supplied mixed with sucrose or mannitol to maintain stability and aid their solubility. For these two proteins, it is particularly important to dissolve the full content of the vial to get an homogeneous solution. It is recommended that carbonic anhydrase is dissolved in distilled water to avoid the formation of aggregates during freezing and thawing.

#### 5) Protein mix preparation

Dilute the proper combination of calibration kit proteins in the buffer. To obtain peaks with similar heights at 280 nm, use the concentrations in Table 5. The concentrations have been calculated with an assumed applied volume of 0.5% of the geometrical column volume ( $V_c$ ). See appendix for a thorough explanation how to prepare a typical protein mix sample.

**Note:** If precipitation of the proteins occurs upon mixing we recommend brief centrifugation to clarify the protein mixture before applying it to the column.

**Table 5.** Suggested protein concentrations for producing peaks of similar height.

Kit	Protein	Protein concentration
LMW	Aprotinin	3 mg/ml
LMW	Ribonuclease A	3 mg/ml
LMW	Carbonic Anhydrase	3 mg/ml
LMW, HMW	Ovalbumin	4 mg/ml
LMW, HMW	Conalbumin	3 mg/ml
HMW	Aldolase*	4 mg/ml
HMW	Ferritin*	0.3 mg/ml
HMW	Thyroglobulin	5 mg/ml

\*Note: These proteins are supplied mixed with sucrose or mannitol to maintain stability.

#### 6) Size of sample volume

Apply Calibration Kit proteins to the column. To get good resolution, the sample size should not exceed 2% of the geometric column volume,  $\rm V_c.$ 

 $(V_c = r^2 \times \pi \times I)$  where r is radius and I is column length).

#### 7) Determination of elution volume ( $V_e$ )

From the UV curve, determine the elution volumes ( $V_e$ ) for the Calibration Kit proteins by measuring the volume of the eluent from the point of injection to the center of the elution peak, see figure 1.

#### 8) Determination of void volume ( $V_o$ )

The elution volume for Blue Dextran 2000 is equal to the column void volume (V<sub>o</sub>). Prepare a fresh solution of Blue Dextran 2000 (1.0 mg/ml) in the buffer. The rate of solubilization of the Blue Dextran 2000 may be increased by heating the buffer to 50°C before adding the Blue Dextran 2000.

Apply a sample to the column (sample size, 0.5% of the geometric column volume) to determine the void volume ( $V_0$ ). The elution of Blue Dextran can be conveniently monitored at wavelengths of 254, 280 or 620 nm.

It is strongly recommended that the Blue Dextran 2000 is run alone, not mixed with the Calibration Kit or sample proteins, as the fraction of Blue Dextran is broad and may overlap the protein peaks. Always calculate the void volume from the first eluted peak from Blue Dextran.



Figure 1. Elution profiles of Calibration Kit proteins on HiLoad 16/60 Superdex 200 pg column. Elution volumes  $(V_e)$  are found at maximum peak height of each respective protein, see for example, carbonic anhydrase.

#### 9) Preparation of calibration curve

Calculate the  $\mathrm{K}_{\mathrm{av}}$  values for the Calibration Kit proteins using the equation:

$$K_{av} = \frac{V_e - V_o}{V_c - V_o}$$

where  $\rm V_{o}$  = column void volume,  $\rm V_{e}$  = elution volume, and  $\rm V_{c}$  = geometric column volume.

Prepare a calibration curve of  $K_{av}$  versus log molecular weight either on semilogarithmic paper or with a calculation program. It should be possible to fit a curve to the data points, see Figures 2 to 10.

#### 10) Molecular weight determination

Apply the unknown sample (volume 0.1 to 2% of V<sub>c</sub>) and determine the elution volume (V<sub>e</sub>) of the compound of interest. Adjust the concentration of the sample taking into consideration that a sample of 0.5% of V<sub>c</sub> will be diluted 5 to 15-fold during the run.

Calculate the corresponding  $\rm K_{av}$  for the component of interest and determine its molecular weight from the calibration curve prepared using the Calibration Kit proteins.

Note that the molecular weight determinations using the molecular weights of glycoproteins, lipoproteins, non-globular proteins, or other polymers may not correlate well to the calibration curves established for globular proteins by the Calibration Kit proteins. For such compounds, useful information can be obtained by relating their elution volume data to a molecular size parameter, such as Stoke's radius ( $R_{st}$ ), rather than to molecular weight values. Plots of  $\sqrt{-\log (K_{av})}$  vs.  $R_{st}$  have been used successfully to determine the Stoke's radius of proteins and our Calibration Kit proteins may be used for these plots too.

### 6. Typical results

#### Method used for figures 2 to 10

Sample:	Proteins from Gel Filtration Calibration Kits LMW and HMW:		
aprotinin (Apr), RNAse A (R), carbonic anhydras ovalbumin (O), conalbumin (C), aldolase (Ald), fe and thyroglobulin (T)			
Sample vol.	: Figures 2 to 5:100 µl Figures 6 to 10:500 µl		
Buffer:	50 mM phosphate buffer, 150 mM NaCl, pH 7.2		
Flow rate:	Figures 2, 4, 5, 8, 9 and 10: Figure 3: Figures 6 and 7:	0.5 ml/min 0.6 ml/min 1.0 ml/min	
System:	ÄKTAexplorer 10		
Detection:	280 nm		
mAU		1.00 Kov	



Figure 2. Chromatographic separation and calibration curve for some of the standard proteins on Superdex 200 10/300 GL column.



Figure 3. Chromatographic separation and calibration curve for some of the standard proteins on Superdex 75 10/300 GL column.



Figure 4. Chromatographic separation and calibration curve for some of the standard proteins on Superose 6 10/300 GL column.

**Note:** Thyroglobulin may be excluded from the calculation of  $K_{av}$  due to none-linear behavior of thyroglobulin on this column. Thyroglobulin may however, be included in a plot of  $\sqrt{-\log (K_{av})}$  vs. Stoke's Radius ( $R_{st}$ ).



**Figure 5.** Chromatographic separation and calibration curve for some of the standard proteins on **Superose 12 10/300 GL** column.



**Figure 6.** Chromatographic separation and calibration curve for some of the standard proteins on **HiLoad 16/60 Superdex 200 pg** column.



**Figure 7.** Chromatographic separation and calibration curve for some of the standard proteins on **HiLoad 16/60 Superdex 75 pg** column.



**Figure 8.** Chromatographic separation and calibration curve for some of the standard proteins on **HiPrep 16/60 Sephacryl S-300 HR** column.

Note: Aprotinin may be excluded from the calculation of  $\rm K_{av}$  due to none-linear behavior of aprotinin on this column.



**Figure 9.** Chromatographic separation and calibration curve for some of the standard proteins on **HiPrep 16/60 Sephacryl S-200 HR** column.



**Figure 10.** Chromatographic separation and calibration curve for some of the standard proteins on **HiPrep 16/60 Sephacryl S-100 HR** column.

### 7. Important notes

# 7.1. The use of the calibration kits with denaturing solvents

The molecular weight ranges given in Table 3 are for globular proteins in their native conformations. The use of denaturing agents, such as sodium dodecyl sulfate (SDS), chaotropic salts and guanidine hydrochloride (GuHCl) and hydrogen bond disrupting agents, such as urea, may alter the molecular conformation of proteins often greatly increasing their hydrodynamic volumes. Since separations by gel filtration are based on molecular size, the molecular weight ranges change when the proteins assume extended conformations.

Superdex 200 has the most useful molecular weight range and flow properties in solvents where proteins are completely denatured (exclusion limit is approximately 120 000 for completely denatured proteins). The Low Molecular Calibration Kit is suitable for the calibration of columns in denaturing solvents. Each protein in the kit comprises of a single polypeptide chain therefore, their molecular weights do not change when they are exposed to denaturants (although their Stoke's Radii do change).

# 7.2. Dimer and oligomer formation in calibration kit proteins

The ribonuclease A, conalbumin, aldolase, ferritin and thyroglobulin standards may contain small amounts of apparent dimers or oligomers that elute in the void volume or slightly before the true peak. Dimers can be used to produce more calibration points, however, pure dimer formation has to be determined before calculations.

### 7.3. Electrophoresis calibration kits

Gel Filtration Calibration Kits (HMW and LMW) contain protein standards for use in gel filtration chromatography only. Kits containing protein standards for molecular weight determination by polyacrylamide gel electrophoresis are also available from GE Healthcare. Please visit, www.gelifesciences.com

## 8. Appendix

#### Protein mix preparation for calibrating HiLoad Superdex 200 pg 26/60

Geometric column volume ( $V_c$ ) = 320 ml Recommended sample volume = 1.6 ml (0.5%  $V_c$ )

#### Preparation of "Mix a" in Table 4 (Ferritin, Conalbumin, Carbonic Anhydrase and Ribonucelase A) for Superdex 200

Prepare each individual protein as described in section 5.4 to a concentration of 20 mg/ml.

The recommended protein concentrations for the different proteins in a mix are according to Table 5:

Ferritin: 0.3 mg/ml Conalbumin: 3 mg/ml Carbonic Anhydrase: 3 mg/ml Ribonuclease A: 3 mg/ml

To prepare a 2 ml protein mix, make the following calculations:

Calculate the volume (X) of each protein needed from the 20 mg/ml solution in order to achieve the desired protein concentration

Ferritin: 20 mg/ml \* X ml = 2 ml \* 0.3 mg/ml X = 0.03 ml

All other proteins: 20 mg/ml \* X ml = 2 ml \* 3 mg/ml X = 0.3 ml

Finally mix 0.03 ml Ferritin, 0.3 ml Conalbumin, 0.3 ml Carbonic Anhydrase and 0.3 ml Ribonucelase and dilute with 1.07 ml buffer to get 2 ml final volume.

Inject 1.6 ml of the protein mix onto the column.

### 9. Product information

Gel Filtration Calibration Kits	Code no.
Low molecular weight	28-4038-41
High molecular weight	28-4038-42

Related products	Quantity	Code no.
Superdex 75 10/300 GL	1	17-5174-01
Superdex 75 5/150 GL	1	28-9205-04
Superdex 75 PC 3.2/30	1	17-0771-01
Superdex 200 10/300 GL	1	17-5175-01
Superdex 200 5/150 GL	1	28-9065-61
Superdex 200 PC 3.2/30	1	17-1089-01
HiLoad 16/60 Superdex 75 pg	1	17-1068-01
HiLoad 26/60 Superdex 75 pg	1	17-1070-01
HiLoad 16/60 Superdex 200 pg	1	17-1069-01
HiLoad 26/60 Superdex 200 pg	1	17-1071-01
Superose 12 10/300 GL	1	17-5173-01
Superose 12 PC 3.2/30	1	17-0674-01
Superose 6 10/300 GL	1	17-5172-01
Superose 6 PC 3.2/30	1	17-0673-01
HiPrep 16/60 Sephacryl S-100 HR	1	17-1165-01
HiPrep 26/60 Sephacryl S-100 HR	1	17-1194-01
HiPrep 16/60 Sephacryl S-200 HR	1	17-1166-01
HiPrep 26/60 Sephacryl S-200 HR	1	17-1195-01
HiPrep 16/60 Sephacryl S-300 HR	1	17-1167-01
HiPrep 26/60 Sephacryl S-300 HR	1	17-1196-01
Precision Column Holder	1	17-1455-01
Reference literature		
Selection guide:		
Gel Filtration Columns and Media		18-1124-19
Handbook:		
Gel Filtration - Principles & Metho	ds	18-1022-18

#### GE Healthcare offices:

GE Healthcare Bio-Sciences AB Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp. 800 Centennial Avenue, P.O.Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

For contact information for your local office, please visit: www.gelifesciences.com/contact GE Healthcare UK Limited Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK http://www.gelifesciences.com/protein-purification



28951560