

## PfuS Mix Purification prototol

Adapted from Görlich laboratory

### I. Expression

	PfuS Polymerise	Pab-PPase	Pab-dTUPase
Overnight culture	pSF421 in BL21DE3Rosetta, 50 ml 2YT with Amp and CAM. 37 °C	pSF318 in BL21DE3Rosetta, 50 ml 2YT with Amp and CAM. 37 °C	pSF323 in BL21DE3Rosetta, 50 ml 2YT with 2% glycerin, Amp and CAM. 37 °C
Dilution	Add 20 ml of overnight culture in 1L prewarmed TB medium. grow until OD=0,5 at 37 °C		Add 20 ml of overnight culture in 1L prewarmed <b>2YT medium with 0.4 M NaCl</b> . grow until OD=0,5 at 37 °C
Induction	add IPTG (final con. <b>0.2 mM</b> ), shake 4 h at 37 °C	add IPTG (final con. <b>0.2 mM</b> ) shake 4 h at 37 °C	add IPTG (final con. <b>1 mM</b> ) shake 4 h at 37 °C
Harvest	add 1mM PMSF, centrifuge at 4k, 20 min The pellet could be freeze at -20 °C.		

### II. Purification

	PfuS Polymerise	Pab-PPase	Pab-dTUPase
Resuspension	thoroughly suspend the cells in HS-buffer with 5 mM DTT, 5 mM Imidazole, 1 mM PMSF and 1 mM Benzamidine. Use 1 ml for 100 OD, e. g. 10 ml for 1000ml culture of OD=1.	thoroughly suspend the cells in <b>LS-buffer</b> with 5 mM DTT, 5 mM Imidazole, 1 mM PMSF and 1 mM BenzamCidine. Use 1 ml for 100 OD, e. g. 10 ml for 1000ml culture of OD=1.	thoroughly suspend the cells in <b>HS-buffer</b> with 5 mM DTT, 5 mM Imidazole, <b>0.1% TX-100</b> , 1 mM PMSF and 1 mM Benzamidine. Use 1 ml for 100 OD, e. g. 10 ml for 1000ml culture of OD=1.

Lysis	lyse cells with the microfluidizer		
Centrifuge	14K, 20 min, 2 times, transfer the clear supernatant to a new vial.		
Histrap column preparation	connect 3x1ml Histrap column tandemly on Äkta. wash with water and equilibrate with <b>HS-buffer</b> .	connect 3x1ml Histrap column tandemly on Äkta. wash with water and equilibrate with <b>LS-buffer</b> .	connect 3x1ml Histrap column tandemly on Äkta. wash with water and equilibrate with <b>HS-buffer</b> .
Column chromatography*	apply the clear supernatant to the column (1 ml/min), collect the flow through.		
	wash with <b>HS-buffer</b> with 5mM DTT and 5mM imidazole for 5 column volume and wash with <b>LS-buffer</b> with 5mM DTT and 5mM imidazole until the A280 does not change anymore	wash with <b>LS-buffer</b> with 5mM DTT and 5mM imidazole until the A280 does not change anymore	wash with <b>HS-buffer</b> with 5mM DTT and 5mM imidazole until the A280 does not change anymore
	elute with <b>LS-buffer</b> with 5 mM DTT and <b>400mM</b> imidazole		elute with <b>HS-buffer</b> with 5 mM DTT and <b>600mM</b> imidazole
	Pool the peak fractions.		
	1OD <sub>280</sub> = 0.73 mg/ml	1OD <sub>280</sub> = 0.70 mg/ml	1OD <sub>280</sub> = 0.75 mg/ml

Note: The Histrap column shows purple color in the presence of DTT.

### III. Dilution.

Dilute and mix to the final concentration with dilution buffer.

PfuS : 100ng/ul  
PPase: 15ng/ul  
dUTPase: 5ng/ul

IV. Buffer:

HS-Buffer	LS-Buffer	Dilution Buffer
2 M NaCl 50 mM Tris-HCl pH 8.0 10 mM MgCl <sub>2</sub>	290 mM NaCl 44 mM Tris-HCl pH 7.5 4.4 mM MgCl <sub>2</sub>	20 mM Tris-HCl pH 7.5 100 mM KCl 1 mM DTT 0.1 mM EDTA 0.5% NP40 0.5% Tween20 100 µg/ml BSA 50% Glycerol